

ROSE PEN GOLD PLATING SOLUTION 10 ML

SKU: 1735



Independent lab tests prove that TWL's B1RL 5N has the world's most extended lasting Rose finish - a longer resistance to oxidation than any other rose plating formula. Fast and easy to use with no fear of dangerous cyanide fumes. Pen plating allows for precise plating over small areas; highlights, setting prongs, individual chain loops, ring sections or bracelet parts. Pen plating is an easy way to get two and three colour designs. Pen-work requires technicians to be very close to the work surface. Make it a priority to protect yourself during this precision work and use our cyanide-free pen formulations.

PRODUCT DESCRIPTION

Independent lab tests prove that TWL's B1RL 5N has the world's most extended lasting Rose finish - a longer resistance to oxidation than any other rose plating formula. Fast and easy to use with no fear of dangerous cyanide fumes. Pen plating allows for precise plating over small areas; highlights, setting prongs, individual chain loops, ring sections or bracelet parts. Pen plating is an easy way to get two and three colour designs. Pen-work requires technicians to be very close to the work surface. Make it a priority to protect yourself during this precision work and use our cyanide-free pen formulations.

[Find out more!](#)

ADDITIONAL INFORMATION



Weight	0.1 kg
Dimensions	2.5 × 2.5 × 5.5 cm
Gold Content	0.1g
Plating Temperature	Room Temperature
Plating Voltage	8-12v
For Plating On	Gold, Silver, Brass, Rhodium, Palladium and Nickel
Other	Requires a plating pen with tip(s) and pen plating rectifier.